




# DATA SHEET

# Clear Image Singles<sup>®</sup> STERILE Ultrasound Scanning Gel

<p><b>DESCRIPTION</b></p> <ul style="list-style-type: none"> <li>Individually packaged, single patient sterile packets of ultrasound coupling media overwrapped in a peelable pouch</li> <li>For diagnostic and therapeutic imaging applications where a sterile gel is desired</li> <li>Provides optimal sound velocity and acoustic impedance</li> <li>Water-based, water-soluble formula; clear, non-staining</li> <li>Contains no added dyes, perfumes, parabens or surfactants</li> <li>Made in USA</li> </ul> 	<p><b>FUNCTION</b></p> <ul style="list-style-type: none"> <li>Couples diagnostic ultrasound between the transducer and the patient</li> <li>Provides increased acoustic transmission and reduced scattering from air bubbles</li> <li>Easy-to-use, single patient packaging helps prevent cross-contamination and helps reduce the risk of nosocomial infection associated with refilling and reusing ultrasound gel bottles</li> </ul>						
<p><b>USE</b></p> <ul style="list-style-type: none"> <li>To open, peel apart the outer pouch to remove the inner gel packet. Then, tear along the perforation and dispense the gel through the built-in dispensing spout</li> <li>After use, dispose of opened packets according to the facility's protocol</li> <li>Recommended for single patient use where sterile gel is desired</li> </ul>	<p><b>TRANSDUCER COMPATIBILITY</b></p> <ul style="list-style-type: none"> <li>Contains no oils, alcohols, silicones, perfumes, dyes or other ingredients that may cause harm to a transducer</li> </ul>						
<p><b>STERILITY</b></p> <ul style="list-style-type: none"> <li>Sterilization method: Gamma irradiation</li> <li>Sterility assurance level (SAL): 10<sup>-6</sup></li> </ul>	<p><b>CONTRAINDICATIONS</b></p> <ul style="list-style-type: none"> <li>Not intended for use as an electromedical conductor</li> </ul>						
<p><b>ORDERING INFORMATION</b></p> <table border="1"> <thead> <tr> <th><i>Description</i></th> <th><i>Part No.</i></th> <th><i>Contents</i></th> </tr> </thead> <tbody> <tr> <td><b>20 Gram</b> sterile single patient packet</td> <td><b>04-048</b></td> <td>48 packets/box; 6 boxes/case</td> </tr> </tbody> </table>	<i>Description</i>	<i>Part No.</i>	<i>Contents</i>	<b>20 Gram</b> sterile single patient packet	<b>04-048</b>	48 packets/box; 6 boxes/case	<p><b>STORAGE</b></p> <ul style="list-style-type: none"> <li>Store away from direct sunlight to preserve viscosity</li> <li>Recommended storage temperature: 32°F to 90°F (0°C to 32°C)</li> <li>Shelf life is three years</li> </ul>
<i>Description</i>	<i>Part No.</i>	<i>Contents</i>					
<b>20 Gram</b> sterile single patient packet	<b>04-048</b>	48 packets/box; 6 boxes/case					