



# DATA SHEET

# SONO<sup>®</sup> Ultrasound Wipes

Distributed by NEXT Medical Products

<p><b>DESCRIPTION</b></p> <ul style="list-style-type: none"> <li>• One-step cleaner, disinfectant, deodorizer and sanitizer specifically formulated for use on ultrasound equipment and transducers.</li> <li>• Thoroughly cleans and disinfects, yet does not damage transducers and equipment.</li> <li>• Kills 99.9% of bacteria and viruses on contact.</li> <li>• Kills 47 pathogens including MRSA, HIV and Staph.</li> <li>• Free of ammonia, phosphate, solvent and bleach; low alcohol level.</li> <li>• Wipes are large enough to provide ample surface area for cleaning ultrasound equipment.</li> <li>• EPA registered. Meets OSHA's Bloodborne Pathogen Standard.</li> <li>• Made in USA</li> </ul>	<p><b>FUNCTION</b></p> <ul style="list-style-type: none"> <li>• One step cleaning, disinfecting and deodorizing wipe for ultrasound equipment.</li> <li>• Can be used on non-porous surfaces as well as transducers, cables, lenses, touch screens and system consoles.</li> <li>• Will not harm most surfaces including acrylic, sealed fiberglass and vinyl.</li> </ul>												
<p><b>USE</b></p> <ul style="list-style-type: none"> <li>• Remove wipe as needed. Consult package for full instructions.</li> <li>• Safe to handle without gloves.</li> <li>• To retain moisture, reseal package firmly.</li> </ul>	<p><b>CONTRAINDICATIONS</b></p> <ul style="list-style-type: none"> <li>• Do not use as a diaper wipe or for personal cleansing. Not a baby wipe.</li> <li>• Do not flush in toilet.</li> <li>• Spot test on plastic or painted surfaces in an inconspicuous area before use.</li> <li>• Do not use on unpainted wood.</li> <li>• Consult package for full product information.</li> <li>• Keep out of reach of children.</li> </ul>												
<p><b>ORDERING INFORMATION</b></p> <table border="1"> <thead> <tr> <th>Description</th> <th>Part No.</th> <th>Wipe Size</th> <th>Contents</th> </tr> </thead> <tbody> <tr> <td>SONO<sup>®</sup> Ultrasound Wipes, soft pack</td> <td>4018</td> <td>7" x 10"</td> <td>50 ultrasound wipes/pack, 12 packs/case</td> </tr> <tr> <td>SONO<sup>®</sup> Ultrasound Wipes, canister</td> <td>4032</td> <td>7" x 8"</td> <td>80 ultrasound wipes/canister, 6 canisters/case</td> </tr> </tbody> </table>	Description	Part No.	Wipe Size	Contents	SONO <sup>®</sup> Ultrasound Wipes, soft pack	4018	7" x 10"	50 ultrasound wipes/pack, 12 packs/case	SONO <sup>®</sup> Ultrasound Wipes, canister	4032	7" x 8"	80 ultrasound wipes/canister, 6 canisters/case	<p><b>STORAGE</b></p> <ul style="list-style-type: none"> <li>• Store in original container in areas inaccessible to children.</li> <li>• Keep package securely closed to retain moisture.</li> </ul>
Description	Part No.	Wipe Size	Contents										
SONO <sup>®</sup> Ultrasound Wipes, soft pack	4018	7" x 10"	50 ultrasound wipes/pack, 12 packs/case										
SONO <sup>®</sup> Ultrasound Wipes, canister	4032	7" x 8"	80 ultrasound wipes/canister, 6 canisters/case										



**Distributed by:**

NEXT Medical Products Company, 45 Columbia Road, Branchburg, NJ 08876  
 Order Desk: 800-458-4254 Order Fax: 908-218-3715 Website: [www.nextmedicalproducts.com](http://www.nextmedicalproducts.com)