

# SONO® Ultrasound Wipes


## SPECIALIZED CARE FOR YOUR ULTRASOUND EQUIPMENT



Most general disinfectant wipes can cause damage to ultrasound equipment, leading to costly repairs or replacement.

SONO Ultrasound Wipes are a proprietary **one-step cleaner, disinfectant, deodorizer and sanitizer** specifically formulated for use on ultrasound equipment and transducers. It thoroughly cleans and disinfects without causing harm to your equipment. Protect your investment with SONO Ultrasound Wipes.

- ✓ Tested on major brands of ultrasound equipment.
- ✓ Kills 47 pathogens including MRSA, HIV and Staph.
- ✓ Safe to handle without gloves.
- ✓ Won't damage equipment; safe for LCD/ touch screens.
- ✓ Helps avoid costly damage to equipment such as delamination of transducer heads and yellowing and cracking of cables.
- ✓ Free of phosphates, ammonia, solvent and bleach; low alcohol level.
- ✓ Convenient - no need for multiple cleaning products. Can be used on non-porous surfaces as well as transducers, cables, lenses, touch screens and system consoles.



*Delamination caused gel to seep inside.*

Using wipes NOT designed for ultrasound equipment can cause serious damage.

Approximate cost to hospital: \$4,500  
These touch screens were cleaned with hydrogen peroxide wipes, causing the panels to become unresponsive.

Replacement Cost \$9,000-\$15,000 can lead to patient shock and image quality issues.

Specialty Transducers can cost upwards of \$25,000  
Alcohol and solvents can dry and crack cables, exposing wires.

**SONO Ultrasound Wipes clean your equipment without harm.**



Two packaging options are available.

Canister:  
P/N 4032  
Wipes are 7" x 8"  
packaged  
80 per  
container, 6/  
case

Soft Pack: P/N 4018  
Wipes are 7" x 10"  
Package of 50, 12  
packs/case



Made *SPECIFICALLY* for Ultrasound!

SONO®

ULTRASOUND WIPES

One Step Cleaner-Disinfectant-Deodorizer-Sanitizer

Distributed by:  
**NEXT Medical Products Company**  
For more information or to place an order, call  
1-800-458-4254 or email [customerservice@nextmedicalproducts.com](mailto:customerservice@nextmedicalproducts.com)

Patent Pending  
Made in the USA

Meets OSHA's Bloodborne Pathogen Standard 1910.1030  
SONO is a registered trademark of Advanced Ultrasound Solutions, Inc.