

METAL MAX™ Reliability Challenge

Please contact NEXT Medical Products Company to take part!
See details on reverse - terms and conditions apply



Combi™ Laryngoscopes

No Compromise Professional Instruments



proactmedical.co.uk

Metal Max™ + Disposable Combi™ Laryngoscope Sets - Preloaded 6V Battery Pack



PROACT Combi™ Sets offer unparalleled reliability and performance and are ready to use in less than 5 seconds.

Now also available in Lever Tip configuration, the Combi™ Laryngoscope range provides an all-in-one solution, including permanently attached handle, preloaded 6V battery pack and complete with a colour coded test-in-bag protective sheath. The simplest, fastest and most dependable solution available - please try for yourself.

We believe that this system is up to 10x more reliable than many other comparable competitive products!

- **Ready for use in less than 5 seconds!**
- **Essential test-in-bag capability for repeated testing and real peace of mind**
- Non-removable Metal Max™ blade permanently attached to a lightweight, super strong, aluminium full metal handle removes hassle and compatibility issues
- Super strong Metal Max™ blade with steel through-the-hook for maximum strength
- Full length blade spatula and ET tube guide surface - no compromise
- 31% lighter than the equivalent set when fitted with normal batteries
- Long battery life - over 5 hours continuous illumination. 3 years battery shelf life from the date of manufacture
- Fast and easy battery removal at end of life
- Eco friendly design ensures ease of disposal and environmental efficiency
- Supplied in easy-to-open colour coded packaging
- New extended range including Stubby handle and Cardiff PRO variants



	Part No.	Fitting	Size	Length (mm)	Length (inch)
	MMCB90MAC0	Combi™	MAC 0	61.00	2.40
	MMCB90MAC1	Combi™	MAC 1	75.00	2.95
	MMCB90MAC2	Combi™	MAC 2	93.00	3.66
	MMCB90MAC3	Combi™	MAC 3	110.00	4.33
	MMCB90MAC4	Combi™	MAC 4	135.00	5.31
	MMCB90MAC5	Combi™	MAC 5	156.00	6.14
	MMCB90MIL00	Combi™	MIL 00	43.50	1.71
	MMCB90MIL0	Combi™	MIL 0	55.00	2.16
	MMCB90MIL1	Combi™	MIL 1	80.50	3.17
	MMCB90MIL2	Combi™	MIL 2	133.00	5.24
	MMCB90MIL3	Combi™	MIL 3	170.00	6.72
	MMCB90MIL4	Combi™	MIL 4	182.00	7.17
	MMCB90LT2	Combi™	Lever Tip 2	93.00	3.66
	MMCB90LT3	Combi™	Lever Tip 3	110.00	4.33
	MMCB90LT4	Combi™	Lever Tip 4	135.00	5.31

	Part No.	Fitting	Size	Length (mm)	Length (inch)
	MMCB90CP00	Combi™ Cardiff PRO	CP 00	55.00	2.17
	MMCB90CP0	Combi™ Cardiff PRO	CP 0	71.00	2.81
	MMCB90CP01	Combi™ Cardiff PRO	CP 1	77.00	3.03
	MMCB90LTCP1	Combi™ Cardiff PRO Lever Tip	CP LT 1	77.00	3.03
	MMCB90MAC3S	Combi™ Stubby Handle	MAC 3	110.00	4.33
	MMCB90MAC4S	Combi™ Stubby Handle	MAC 4	135.00	5.31

(Available late 2017)



Lever Tip Combi™ Mac Blade

All blades size 2 and smaller are fitted with a paediatric slimline handle with battery

ESSENTIAL TEST IN BAG PACKAGING



North and South American Distribution by:
 NEXT Medical Products Company
 45 Columbia Road, Branchburg, NJ 08876, USA
 Tel: 908-722-4549 or 800-458-4254 Fax: 908-218-3715
 e: customerservice@nextmedicalproducts.com
 website: www.nextmedicalproducts.com



Manufactured by:

PROACT Medical Ltd
 9-13 Oakley Hay Lodge
 Great Folds Road
 Oakley Hay Business Park
 Corby, Northants, NN18 9AS
 United Kingdom

Tel: 01536 461981
 Fax: 01536 461837
 Intl Tel: +44(0) 1536 461981
 Intl Fax: +44(0) 1536 461837
 e: sales@proactmedical.co.uk
 www.proactmedical.co.uk